

m-PEG3-aminoxy-Boc

Chemical Properties

CAS No. : 1835759-87-9


Formula: C₁₂H₂₅N₁O₆

Molecular Weight: 279.33

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-PEG3-ONHBoc is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands linked by a connector: one ligand binds to an E3 ubiquitin ligase, while the other targets a specific protein. These compounds leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.580 mL	17.900 mL	35.7999 mL
5 mM	0.716 mL	3.580 mL	7.160 mL
10 mM	0.358 mL	1.790 mL	3.580 mL
50 mM	0.0716 mL	0.358 mL	0.716 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481