


## Licoricone

## Chemical Properties

CAS No. : 51847-92-8  
 Formula: C<sub>22</sub>H<sub>22</sub>O<sub>6</sub>  
 Molecular Weight: 382.41  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Licoricone exhibits anti-helicobacter pylori activity against the CLAR and AMOX-resistant strain as well as four CLAR (AMOX)-sensitive strains.
Targets(IC50)	Antibacterial,Antifection

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.615 mL	13.075 mL	26.1499 mL
5 mM	0.523 mL	2.615 mL	5.230 mL
10 mM	0.2615 mL	1.3075 mL	2.615 mL
50 mM	0.0523 mL	0.2615 mL	0.523 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481