


## Ydspstst peptide

## Chemical Properties

CAS No. :	130007-45-3
Formula:	C <sub>38</sub> H <sub>57</sub> N <sub>9</sub> O <sub>19</sub>
Molecular Weight:	943.91
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Ydspstst peptide is a biochemical.
-------------	------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0594 mL	5.2971 mL	10.5942 mL
5 mM	0.2119 mL	1.0594 mL	2.1188 mL
10 mM	0.1059 mL	0.5297 mL	1.0594 mL
50 mM	0.0212 mL	0.1059 mL	0.2119 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481