

ICI-147798 HCl

## Chemical Properties


CAS No. : 86015-39-6

Formula: C<sub>24</sub>H<sub>31</sub>Cl<sub>2</sub>N<sub>7</sub>O<sub>4</sub>

Molecular Weight: 552.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	ICI-147798 HCl is a beta adrenoceptor antagonist.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8101 mL	9.0506 mL	18.1012 mL
5 mM	0.362 mL	1.8101 mL	3.6202 mL
10 mM	0.181 mL	0.9051 mL	1.8101 mL
50 mM	0.0362 mL	0.181 mL	0.362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481