

Compound T127792

Chemical Properties

CAS No. :

Formula: C17H30O14

Molecular Weight: 458.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2S,3R,4R,5R,6R)-2-methyl-6-(((2R,3S,4S,5R,6S)-3,4,6-trihydroxy-5-(((2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxy)oxan-2-yl]methoxy)oxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1815 mL	10.9073 mL	21.8145 mL
5 mM	0.4363 mL	2.1815 mL	4.3629 mL
10 mM	0.2181 mL	1.0907 mL	2.1815 mL
50 mM	0.0436 mL	0.2181 mL	0.4363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481