


L803-mts

Chemical Properties

CAS No. : 1043881-55-5
 Formula: C₆₆H₁₁₀N₁₅O₂₀P
 Molecular Weight: 1464.665
 Storage: Keep away from moisture
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L803-mts has a wide range of applications in life science related research.
Targets(IC50)	Beta Amyloid,GSK-3

Solubility Information

Solubility	20% Acetonitrile: 1 mg/mL (0.68 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6827 mL	3.4137 mL	6.8275 mL
5 mM	0.1365 mL	0.6827 mL	1.3655 mL
10 mM	0.0683 mL	0.3414 mL	0.6827 mL
50 mM	0.0137 mL	0.0683 mL	0.1365 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481