

Compound T130293

Chemical Properties


CAS No. :

Formula: C54H90O30

Molecular Weight: 1219.28

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-({6-[(6-[(6-[(2-[(3,4-dihydroxy-5-[(3,4,5-trihydroxyoxan-2-yl)oxy]oxan-2-yl)oxy]methyl]-3,5-dihydroxy-6-[(2E,6E)-3,7,11-trimethyldodeca-2,6,10-trien-1-yl]oxy)oxan-4-yl]oxy]-4,5-dihydroxy-2-methyloxan-3-yl]oxy]-4,5-dihydroxy-2-methyloxan-3-yl]oxy]-4,5-dihydroxyoxan-3-yl]oxy)-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8202 mL	4.1008 mL	8.2016 mL
5 mM	0.164 mL	0.8202 mL	1.6403 mL
10 mM	0.082 mL	0.4101 mL	0.8202 mL
50 mM	0.0164 mL	0.082 mL	0.164 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481