


RMI-12330A

Chemical Properties

CAS No. : 82985-31-7
 Formula: C₂₃H₃₆N₂
 Molecular Weight: 340.55
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	RMI-12330A is an adenylyl cyclase inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9364 mL	14.6821 mL	29.3643 mL
5 mM	0.5873 mL	2.9364 mL	5.8729 mL
10 mM	0.2936 mL	1.4682 mL	2.9364 mL
50 mM	0.0587 mL	0.2936 mL	0.5873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481