


γ -Secretase modulator 4

Chemical Properties

CAS No. :	1420200-82-3
Formula:	C ₂₃ H ₁₉ FN ₄ O ₂
Molecular Weight:	402.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	γ -Secretase modulator 4 is a potent γ -secretase modulator that reduces A β 42 levels with IC ₅₀ s of 0.017 μ M in mice and 0.014 μ M in humans.
Targets(IC ₅₀)	Others,Gamma-secretase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.485 mL	12.4248 mL	24.8497 mL
5 mM	0.497 mL	2.485 mL	4.9699 mL
10 mM	0.2485 mL	1.2425 mL	2.485 mL
50 mM	0.0497 mL	0.2485 mL	0.497 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Bischoff FP, et al. Design and synthesis of a novel series of cyanoindole derivatives as potent γ -secretase modulators. *Bioorg Med Chem Lett.* 2019 Jul 15;29(14):1737-1745.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481