


Ethyl 4-aminophenylacetate

Chemical Properties

CAS No. :	5438-70-0
Formula:	C ₁₀ H ₁₃ NO ₂
Molecular Weight:	179.219
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl 4-aminophenylacetate, with CAS No. 5438-70-0, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Ethyl 4-aminophenylacetate provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.5797 mL	27.8987 mL	55.7973 mL
5 mM	1.1159 mL	5.5797 mL	11.1595 mL
10 mM	0.558 mL	2.7899 mL	5.5797 mL
50 mM	0.1116 mL	0.558 mL	1.1159 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481