

Fmoc-Gln(Trt)-Thr(psi(Me,Me)pro)-OH

Chemical Properties

CAS No. : 1572725-72-4


Formula: C₄₆H₄₅N₃O₇

Molecular Weight: 751.87

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fmoc-Gln(Trt)-Thr(psi(Me,Me)pro)-OH is a dipeptide.
Targets(IC50)	Amino Acids and Derivatives

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.330 mL	6.6501 mL	13.3002 mL
5 mM	0.266 mL	1.330 mL	2.660 mL
10 mM	0.133 mL	0.665 mL	1.330 mL
50 mM	0.0266 mL	0.133 mL	0.266 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481