

## TM4SF2/TSPAN7 Protein, Rat, Recombinant (His)

### General Information

Synonyms:	tetraspanin 7
Protein Construction:	A DNA sequence encoding the rat TSPAN7 (ABX10436.1) (Arg113-Met213) was expressed with a polyhistidine tag at the N-terminus. Predicted N terminal: His
Species:	Rat
Expression Host:	HEK293 Cells
Accession:	B3STU1
Molecular Weight:	13.9 kDa (predicted)

### QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 $\mu$ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

**Reconstitution:**  
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

**Stability & Storage:**

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

**Shipping:**

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

TALLA-1, also known as TSPAN7, is a member of the transmembrane 4 superfamily. Most members of this family are cell-surface proteins that are characterized by the presence of four hydrophobic domains. TALLA-1 gene is associated with X-linked mental retardation and neuropsychiatric diseases such as Huntington's chorea, fragile X syndrome and myotonic dystrophy. TALLA-1 is a cell surface glycoprotein and may have a role in the control of neurite outgrowth. It is known to complex with integrins.

Reference

Berditchevski F. 2002, J Cell Sci. 114 (23): 4143-51.

Castellví-Bel S. et al., 2001, Mol Genet Metab. 72 (2): 104-8.

Abidi FE. et al., 2002, J Med Genet. 39 (6): 430-3.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481