

B7-H6 Protein, Rat, Recombinant (hFc)

General Information

Synonyms:	B7H6;natural killer cell cytotoxicity receptor 3 ligand 1
Protein Construction:	A DNA sequence encoding the rat NCR3LG1 (XP_006223356.1) (Met1-Ser308) was expressed with the Fc region of human IgG1 at the C-terminus. Predicted N terminal: Met
Species:	Rat
Expression Host:	HEK293 Cells
Accession:	XP_006223356.1
Molecular Weight:	60.4 kDa (predicted)

QC Testing

Biological Activity:	1.Measured by its binding ability in a functional ELISA. 2. Immobilized rat NCR3-His at 10 µg/mL (100µL/well) can bind rat B7-H6-Fc , the EC50 of rat B7-H6 -Fc is 1-20ng/mL.
Purity:	> 95 % as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481