

## CD5L Protein, Mouse, Recombinant (His)

### General Information

|                       |  |
|-----------------------|--|
| Synonyms:             | 42010;AAC-11;Sp-alpha;Pdp;Api6;AI047839;CT2;CD5 molecule-like;Sp- $\alpha$ |
| Protein Construction: | Glu22-Val352   |
| Species:              | Mouse  |
| Expression Host:      | HEK293 Cells   |
| Accession:            | Q9QWK4   |
| Molecular Weight:     | 37.7 kDa (Predicted); 45-60 kDa (Due to glycosylation)                     |

### QC Testing

|                      |  |
|----------------------|--|
| Biological Activity: | Activity testing is not tested. It is theoretically active, but we cannot guarantee it.  |
| Purity:              | > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC  |
| Endotoxin:           | < 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.  |
| Formulation:         | Lyophilized from 0.22 $\mu$ m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization. |

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

CD5L, a soluble protein belonging to the SRCR superfamily, is expressed mostly by macrophages in lymphoid and inflamed tissues. The expression of this protein is transcriptionally controlled by LXRs, members of the nuclear receptor family that play major roles in lipid homeostasis. Research undertaken over the last decade has uncovered critical roles of CD5L as a PRR of bacterial and fungal components and in the control of key mechanisms in inflammatory responses, with involvement in processes, such as infection, atherosclerosis, and cancer.

Reference

Resnick, D. et al., 1994, Trends Biochem. Sci. 19: 5-8.

Gebe, J. A. et al., 1997, J. Biol. Chem. 272 (10): 6151-6158.

Sarrias, M.R. et al., 2004, Crit. Rev. Immunol. 24: 1-37.

Mukhopadhyay, S. and Gordon, S., 2004, Immunobiology 209: 39-49.

Sarrias, M. R. et al., 2005, J. Biol. Chem. 280 (42): 35391-35398.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481