

ULBP-1 Protein, Human, Recombinant (His)

General Information

Synonyms:	UL16 binding protein 1;RAET11
Protein Construction:	A DNA sequence encoding the mature form of human ULBP1 (NP_079494.1) (Met 1-Gly 216) was expressed, with a polyhistidine tag at the C-terminus. Predicted N terminal: Gly 26
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q9BZM6
Molecular Weight:	23.8 kDa (predicted); 28-32 kDa (reducing condition, due to glycosylation)

QC Testing

Biological Activity:	Immobilized human ULBP1-His at 10 µg/ml (100 µl/well) can bind human S4-Fc3L3-NKG2D , The EC50 of human S4-Fc3L3-NKG2D is 0.39-0.91 µg/ml.
Purity:	> 98 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

UL16-binding proteins (ULBP) or retinoic acid early transcripts-1 (RAET1) are ligands to the activating receptor, NKG2D. Ten members of the human ULBP/RAET1 gene family have been identified to encode for potentially functional proteins, and have tissue-specific expressions. ULBP1, also known as RAET1I and NKG2DL1, together with at least ULBP 2 and 3, are well-known ligands for NKG2D, and activate multiple signaling pathways in primary NK cells, resulting in the production of cytokines and chemokines. ULBP1 is expressed in T-cells, B-cells,

erythroleukemia cell lines and in a wide range of tissues including heart, brain, lung, liver and bone marrow, as well as some tumor cells. As an unconventional member of the MHC class I family, ULBP1 function in immune responses, especially in cancer and infectious diseases. Unlike other ULBP members, ULBP1 is able to interact with soluble CMV glycoprotein UL16 in CMV infected cells. The interaction with UL16 blocked the interaction with the NKG2D receptor, and thus might escape the immune surveillance. Furthermore, UL16 also causes ULBP1 to be retained in the ER and cis-Golgi apparatus so that it does not reach the cell surface. The ULBP1 regulation may have implications for development of new therapeutic strategies against cancer cells.

Reference

- Rölle, A. et al., 2003, J Immunol. 171(2): 902-908.2.
López-Soto, A. et al., 2006, J Biol Chem. 281(41): 30419-30430.3.
Song, H. et al., 2006, Cell Immunol. 239(1): 22-30.4.
Eisele, G. et al., 2006, Brain. 129 (9): 2416-2425.5.
Romphruk, AV. et al., 2009, Immunogenetics. 61(9): 611-617.6.
Sutherland, C.L. et. al., 2002, J. Immunol. 168: 671-679.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481