

NOV/CCN3 Protein, Human, Recombinant (His)

General Information

Synonyms:	IGFBP-9;CCN3;IGFBP9;NOVh;IBP-9;nephroblastoma overexpressed
Protein Construction:	A DNA sequence encoding the full length of human NOV (NP_002505.1) (Met 1-Met 357) was expressed, fused with a polyhistidine tag at the C-terminus. Predicted N terminal: Thr 32
Species:	Human
Expression Host:	Baculovirus Insect Cells
Accession:	P48745
Molecular Weight:	37 kDa (predicted); 47 kDa (reducing conditions)

QC Testing

Biological Activity:	Measured by the ability of the immobilized protein to support the adhesion of Balb/3T3 mouse embryonic fibroblast cells. When cells are added to CCN3-coated plates (10 µg/ml, 100 µl/well), >50% cells will adhere specifically after 60 minutes at 37 °C.
Purity:	> 94 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 50 mM Tris, 100 mM NaCl, 0.5 mM PMSF, 10 mM Imidazole, 10% Glycerol, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Protein NOV homolog, also known as Nephroblastoma-overexpressed gene protein homolog, NOV, and CCN3, is a putative ligand for integrin receptors, is tightly associated with the extracellular matrix, and mediates diverse cellular functions, including cell adhesion and proliferation. CCN3 has been shown to negatively regulate growth

although it promotes migration in a cell type-specific manner. This secreted protein belongs to the CCN family, and its expression was observed in a broad variety of tissues from the early stage of development, and altered expression of CCN3 has been observed in a variety of tumors, including hepatocellular carcinomas, Wilm's tumors, Ewing's sarcomas, gliomas, rhabdomyosarcomas, and adrenocortical carcinomas. Mature CCN3 protein has five distinct modules and truncated protein variants with altered function are found in many cancers. CCN3 acts through the core stem cell signaling pathways including Notch and Bone Morphogenic Protein, connecting CCN3 with the modulation of self-renewal and maturation of some cell lineages including hematopoietic, osteogenic, and chondrogenic. CCN3 may affect the extracellular environment of the niche for hematopoietic stem cells. CCN3 has emerged as a key player in stem cell regulation, hematopoiesis, and a crucial component within the bone marrow microenvironment.

Reference

- Manara MC, et al. (2002) The expression of ccn3(nov) gene in musculoskeletal tumors. *Am J Pathol.* 160(3): 849-59.
- Lin CG, et al. (2003) CCN3 (NOV) is a novel angiogenic regulator of the CCN protein family. *J Biol Chem.* 278(26): 24200-8.
- Vallacchi V, et al. (2009) CCN3/nephroblastoma overexpressed matricellular protein regulates integrin expression, adhesion, and dissemination in melanoma. *Cancer Res.* 68(3): 715-23.
- Sin WC, et al. (2009) Matricellular protein CCN3 (NOV) regulates actin cytoskeleton reorganization. *J Biol Chem.* 284(43): 29935-44.
- McCallum L, et al. (2009) CCN3--a key regulator of the hematopoietic compartment. *Blood Rev.* 23(2): 79-85.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481