

## Influenza A H1N1 (A/Puerto Rico/8/34/Mount Sinai) Nucleoprotein/NP Protein (I116M, His)

### General Information

Synonyms:	NP Protein
Protein Construction:	A DNA sequence encoding the Influenza A virus (A/Puerto Rico/8/34/Mount Sinai (H1N1)) nucleoprotein (AAM75159.1, with mutation Ile 116 Met) (Met1-Gly490) was fused with a polyhistidine tag at the C-terminus. Predicted N terminal: Met 1
Species:	H1N1
Expression Host:	Baculovirus Insect Cells
Accession:	P03466
Molecular Weight:	56.7 kDa (predicted); 50 kDa (reducing conditions)

### QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 60 mM Tris, 500 mM NaCl, 10% glycerol, pH 7.5. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

Reconstitution:	Reconstituted with sterile deionized water to 0.25 mg/mL. Reconstitution conditions may vary depending on the lot.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>

### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Influenza viral nucleoprotein (NP) plays a critical role in virus replication and host adaptation. Nucleoprotein (NP) is an essential multifunctional protein that encapsidates the viral genome and functions as an adapter between the virus and the host cell machinery. NPs contain two nuclear localization signals (NLSs): a well-studied monopartite NLS1 and a less-characterized NLS2, thought to be bipartite. The nucleocapsid is a complex of the

## A DRUG SCREENING EXPERT

---

viral nucleoprotein, RNA, and several other viral proteins. The nucleoprotein forms large, RNA-bound, helical filaments and acts as a scaffold for additional viral proteins.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481