

## PVRL1/NECTIN1 Protein, Human, Recombinant (His)

### General Information

Synonyms:	HVEC;ED4;PRR1;PVRR1;nectin-1;HIgR;poliovirus receptor-related 1 (herpesvirus entry mediator C);CLPED1;OFC7;PRR;CD111;HV1S;PVRR;SK-12
Protein Construction:	A DNA sequence encoding the human PVRL1 isoform 1 (NP_002846.3) extracellular domain (Met 1-Thr 334) was expressed, with a polyhistidine tag at the C-terminus. Predicted N terminal: Gln 31
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q15223-1
Molecular Weight:	35.4 kDa (predicted); 45-50 kDa (reducing condition, due to glycosylation)

### QC Testing

Biological Activity:	<ol style="list-style-type: none"><li>1. Measured by its binding ability in a functional ELISA. Immobilized PVRL3 at 1 µg/ml (100 µl/well) can bind biotinylated recombinant human PVRL1 / Nectin-1 with a linear range of 6.4-800 ng/ml.</li><li>2. Immobilized Nectin 3 Protein, Human, Recombinant (ECD, hFc Tag) at 2 µg/ml (100 µl/well) can bind PVRL1/NECTIN1 Protein, Human, Recombinant (His Tag), the EC50 of PVRL1/NECTIN1 Protein, Human, Recombinant (His Tag) is 150-700 ng/mL.</li></ol>
Purity:	> 98 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

Reconstitution:	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Stability & Storage:	<p>It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.</p> <p><i>Actual storage temperature shall be subject to the COA.</i></p>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Poliovirus receptor-related 1 (herpesvirus entry mediator C; nectin-1; CD111), also known as PVRL1 is a cell adhesion molecule belonging to the immunoglobulin superfamily that can bind to virion glycoprotein D (gD) to mediate entry of herpes simplex viruses (HSV) and pseudorabies virus (PRV). CD111/Nectin-1/PVRL1 colocalizes with E-cadherin at adherens junctions in epithelial cells. The disruption of cell junctions can result in the redistribution of nectin-1. To determine whether disruption of junctions by calcium depletion influenced the susceptibility of epithelial cells to viral entry, Madin-Darby canine kidney cells expressing endogenous nectin-1 or transfected human nectin-1 were tested for the ability to bind soluble forms of viral gD and to be infected by HSV and PRV, before and after calcium depletion. It has been revealed that binding of HSV and PRV gD was localized to adherens junctions in cells maintained in normal medium but was distributed, along with nectin-1, over the entire cell surface after calcium depletion. Both the binding of gD and the fraction of cells that could be infected by HSV-1 and PRV were enhanced by calcium depletion. Taken together, CD111/Nectin-1/PVRL1 confined to adherens junctions in epithelial cells is not very accessible to virus, whereas dissociation of cell junctions releases nectin-1 to serve more efficiently as an entry receptor.

### Reference

- Yoon M, et al. (2002) Disruption of adherens junctions liberates nectin-1 to serve as receptor for herpes simplex virus and pseudorabies virus entry. *J Virol.* 76(14): 7203-8.
- Mata M, et al. (2001) HveC (nectin-1) is expressed at high levels in sensory neurons, but not in motor neurons, of the rat peripheral nervous system. *J Neurovirol.* 7(5): 476-80.
- HAarr L, et al. (2001) Transcription from the gene encoding the herpesvirus entry receptor nectin-1 (HveC) in nervous tissue of adult mouse. *Virology.* 287(2): 301-9.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481