

Leukotriene A4 Hydrolase Protein, Mouse, Recombinant (His)

General Information

Synonyms:	Lta4h; leukotriene A4 hydrolase
Protein Construction:	A DNA sequence encoding the mouse LTA4H (NP_032543.2) (Met 1-Asp 611) was expressed, with a C-terminal polyhistidine tag. Predicted N terminal: Met 1
Species:	Mouse
Expression Host:	Baculovirus Insect Cells
Accession:	P24527
Molecular Weight:	70.4 kDa (predicted); 62 kDa (reducing conditions)

QC Testing

Biological Activity:	Measured by its ability to cleave the fluorogenic peptide substrate, Arg-7-amido-4-methylcoumarin (R-AMC). The specific activity is >15 pmoles/min/ μ g.
Purity:	> 92 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing 50 mM Tris, 100 mM NaCl, pH 8.0. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Leukotriene A-4 hydrolase, also known as LTA-4 hydrolase, Leukotriene A (4) hydrolase, LTA4H, and LTA4, is a cytoplasm protein that belongs to the peptidase M1 family. LTA4H hydrolyzes an epoxide moiety of leukotriene A4 (LTA-4) to form leukotriene B4 (LTB-4). This enzyme also has some peptidase activity. The leukotrienes (LTs) are a class of structurally related lipid mediators involved in the development and maintenance of inflammatory and allergic reactions. In the biosynthesis of LTs, arachidonic acid was converted into the unstable intermediate

epoxide LTA₄, which may, in turn, be conjugated with glutathione to form the spasmogenic LTC₄, or hydrolyzed into the pro-inflammatory lipid mediator LTB₄ in a reaction catalyzed by Leukotriene A₄ hydrolase (LTA₄H). LTB₄ is a classical chemoattractant of human neutrophils and triggers adherence and aggregation of leukocytes to vascular endothelium, and also modulates immune responses. As a bifunctional zinc metalloenzyme, LTA₄H also exhibits an anion-dependant arginyl aminopeptidase activity of high efficiency and specificity in addition to its epoxide hydrolase activity. LTA₄H is regarded as a therapeutic target for inflammation.

Reference

- Mancini, JA. et al., 1995, Eur. J. Biochem. 231: 65-71.
Orning, L. et al., 1994, J. Biol. Chem. 269: 11269-73.
Rudberg, PC. et al., 2004, J. Biol. Chem. 279: 27376-82.
Qiu, H. et al., 2006, Proc. Natl. Acad. Sci. USA. 103: 8161-6.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481