

Carboxylesterase 2/CES2 Protein, Human, Recombinant (His)

General Information

Synonyms:	carboxylesterase 2;CES2A1;CE-2;PCE-2;iCE
Protein Construction:	A DNA sequence encoding the human CES2 isoform 1 (O00748-1) (Met 1-Leu 559) was expressed, with a C-terminal polyhistidine tag. Predicted N terminal: Gln 27
Species:	Human
Expression Host:	HEK293 Cells
Accession:	O00748-1
Molecular Weight:	60.4 kDa (predicted); 60.4 kDa (reducing conditions)

QC Testing

Biological Activity:	Measured by its ability to hydrolyze p-nitrophenylacetate. The specific activity is >20000 pmols/min/ug.
Purity:	> 95 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 50 mM NaAc, 150 mM NaCl, 10% Glycerol, pH 5. 5. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:	Reconstituted with sterile deionized water to 0.25 mg/mL. Reconstitution conditions may vary depending on the lot.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Carboxylesterase 2 (CES2) is a member of the carboxylesterase family and belongs to the multigene family. Carboxylesterase 2 is responsible for the hydrolysis of ester- and amide-bond-containing drugs such as cocaine and beroin. It also serves to hydrolyze long-chain fatty acid esters and thioesters. It is speculated that carboxylesterases may play a role in lipid metabolism and the blood-brain barrier system and together with isoform 1, are a serine esterase involved in both drug metabolism and activation. Human carboxylesterase 2 is

commonly expressed in tumor tissues and irinotecan, a topoisomerase I inhibitor commonly used in the treatment of many solid tumors.

Reference

Imai T, et al. (1997) Identification and molecular characterization of fractalkine receptor CX3CR1, which mediates both leukocyte migration and adhesion. *Cell*. 91(4): 521-30.

Papadopoulos EJ, et al. (1999) Fractalkine, a CX3C chemokine, is expressed by dendritic cells and is up-regulated upon dendritic cell maturation. *Eur J Immunol*. 29 (8): 2551-9.

Umehara H, et al. (2004) Fractalkine in vascular biology: from basic research to clinical disease". *Arterioscler. Thromb Vasc Biol*. 24 (1): 34-40.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481