

CD32A Protein, Human, Recombinant (167 His, His)

General Information

Synonyms:	IGFR2;CDw32;CD32;CD32A;FCGR2;Fc gamma RIIA;FcGR;Fc γ RIIA;FCG2;FCGR2A1;Fc fragment of IgG receptor IIa
Protein Construction:	A DNA sequence encoding the extracellular domain (Met1-Ile218) of human FCGR2A (P12318-1) was expressed with a C-terminal polyhistidine tag. Predicted N terminal: Gln 34
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P12318-1
Molecular Weight:	22 kDa (predicted); 30 kDa (reducing condition, due to glycosylation)

QC Testing

Biological Activity:	1.Loaded Rituximab (IgG1) on ProA Biosensor, can bind Human Fc γ RIIA / CD32a (H167) recombinant protein with an affinity constant of 1.7 μ M as determined in a BLI assay. 2.Loaded Fc γ RIIA / CD32a (H167) recombinant protein on His1K Biosensor, can bind Bevacizumab (IgG1) with an affinity constant of 0.55 μ M as determined in a BLI assay. 3.Captured Bevacizumab (IgG1) on proA Chip can bind Human Fc γ RIIA / CD32a (H167) recombinant protein with an affinity constant of 0.58 μ M as determined in an SPR assay.
Purity:	> 95 % as determined by SDS-PAGE. > 90 % as determined by SEC-HPLC
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481