

## CD32A Protein, Human, Recombinant (His &amp; Avi), Biotinylated

## General Information

Synonyms:	FCG2;Fc gamma RIIA;FcGR;CD32A;CD32;Fc fragment of IgG receptor IIa;Fc $\gamma$ RIIA;FCGR2;CDw32;IGFR2;FCGR2A1
Protein Construction:	A DNA sequence encoding the extracellular domain (Met1-Met210) of human CD32a (P12318-1) was fused with a c-terminal polyhistidine tagged AVI tag at the C-terminus. The expressed protein was biotinylated in vivo by the Biotin-Protein ligase (BirA enzyme) which is co-expressed. Predicted N terminal: Gln 34
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P12318-1
Molecular Weight:	23.7 kDa (predicted); 33 kDa (reducing condition, due to glycosylation)

## QC Testing

Biological Activity:	1.Loaded Biotinylated Human Fc $\gamma$ RIIA / CD32a (H167) recombinant protein (His & Avi Tag) on SA Biosensor, can bind Rituximab (IgG1) with an affinity constant of 0.53 $\mu$ M as determined in a BLI assay.2. The purity of Human CD32a / FCGR2A Protein (167 His, His & AVI Tag), Biotinylated is more than 95% and the molecular weight of this protein is around 20-38 kDa verified by SEC-MALS.
Purity:	$\geq$ 95 % as determined by SDS-PAGE. $\geq$ 95 % as determined by SEC-HPLC
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 $\mu$ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

## Preparation and Storage

Reconstitution:	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <i>Actual storage temperature shall be subject to the COA.</i>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481