

## CD46 Protein, Human, Recombinant (His)

## General Information

Synonyms:	AHUS2;MIC10;CD46 molecule, complement regulatory protein;TLX;MCP;TRA2.10
Protein Construction:	A DNA sequence encoding the human CD46 (EAW93471.1) extracellular domain (Met 1-Asp 328) was fused with a polyhistidine tag at the C-terminus. Predicted N terminal: Cys 35
Species:	Human
Expression Host:	HEK293 Cells
Accession:	EAW93471.1
Molecular Weight:	34 kDa (predicted); 55-60 kDa (reducing condition, due to glycosylation)

## QC Testing

Biological Activity:	1.Loaded YS5FL derived anti-CD46 antibody, human IgG1 on proA Biosensor, can bind Recombinant Human CD46 Protein, His Tag (Cat#TMPY-02023) with an affinity constant of 0.737 $\mu$ M as determined in BLI assay (Sartorius Octet RED384) (Routinely tested). 2. Immobilized Recombinant Human CD46 Protein (His Tag) (Cat#TMPY-02023) at 1 $\mu$ g/mL (100 $\mu$ L/well) can bind YS5FL derived Anti-CD46 antibody, the EC50 is 3.5-11.0 ng/mL (QC tested).
Purity:	$\geq$ 95 % as determined by SDS-PAGE. $\geq$ 95 % as determined by SEC-HPLC.
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 $\mu$ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

## Preparation and Storage

**Reconstitution:**  
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

**Stability & Storage:**

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

**Shipping:**

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

## Protein Background

CD46, also known as Membrane Cofactor Protein (MCP), is a complement regulatory protein. CD46 is a type 1 membrane protein that plays an important inhibitory role in the complement system. CD46 is expressed in white

blood cells, platelets, epithelial cells, and fibroblasts. Human CD46 shares 5% amino acid sequence identity with mouse and rat CD46. The importance of CD46 to complement regulation is underscored by the observation that genetic loss of CD46 leads to development of atypical hemolytic-uremic syndrome (aHUS), a disease characterized by uncontrolled complement activation. CD46 is implicated in the development and/or progression of selected cancer types.

### Reference

Lublin D.M., et al., (1988), Molecular cloning and chromosomal localization of human membrane cofactor protein (MCP). Evidence for inclusion in the multigene family of complement-regulatory proteins. *J. Exp. Med.* 168:181-194.

Purcell D.F., et al., (1991), Alternatively spliced RNAs encode several isoforms of CD46 (MCP), a regulator of complement activation. *Immunogenetics* 33:335-344.

Post T.W., et al., (1991), Membrane cofactor protein of the complement system: alternative splicing of serine/threonine/proline-rich exons and cytoplasmic tails produces multiple isoforms that correlate with protein phenotype. *J. Exp. Med.* 174:93-102.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 34 Washington Street, Wellesley Hills, MA 02481