

## KIR2DL3 Protein, Human, Recombinant (His)

### General Information

Synonyms:	NKAT2;MGC129943;KIR-K7c;NKAT2B;CU464054.1;NKAT;KIR-023GB;p58;KIR-K7b;KIR2DS5;killer cell immunoglobulin like receptor, two Ig domains and long cytoplasmic tail 3;NKAT2A;GL183;KIRCL23;CD158B2;CD158b
Protein Construction:	A DNA sequence encoding the human KIR2DL3 (AAX23102.1) extracellular domain (Met 1-His 245) was fused with a polyhistidine tag at the C-terminus. Predicted N terminal: His 22
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q32WE4
Molecular Weight:	25.9 kDa (predicted); 45 kDa (reducing condition, due to glycosylation)

### QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 98 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

#### Reconstitution:

Reconstituted with sterile deionized water to 0.25 mg/mL. Reconstitution conditions may vary depending on the lot.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Killer cell immunoglobulin-like receptor 2DL3, also known as CD158 antigen-like family member B2, KIR-23GB, Killer inhibitory receptor cl 2-3, MHC class I NK cell receptor, NKAT2a, NKAT2b, Natural killer-associated transcript 2, p58 natural killer cell receptor clone CL-6, p58.2 MHC class-I-specific NK receptor, CD158b2, and KIR2DL3, is a

single-pass type I membrane protein which belongs to the immunoglobulin superfamily. KIR2DL3 contains 2 Ig-like C2-type (immunoglobulin-like) domains. KIR2DL3 interacts with ARRB2. KIR2DL3 is a receptor on natural killer (NK) cells for HLA-C alleles (HLA-Cw1, HLA-Cw3, and HLA-Cw7). KIR2DL3 inhibits the activity of NK cells thus preventing cell lysis.

### Reference

Selvakumar A., et al., 1997, Immunol. Rev. 155:183-196.

Wilson M.J., et al., 1997, Tissue Antigens 49:574-579.

Maenaka K., et al., 1999, Structure 7:391-398.

Vitale, M. et al., 2004, Int Immunol. 16 (10):1459-66.

Yu M.-C., et al., 2008, Nat. Immunol. 9:898-907.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481