

Carbonic Anhydrase 8 Protein, Mouse, Recombinant (His)

General Information

Synonyms:	RP23-180H12.1;Ca8;AW546993;Cals;Carp;carbonic anhydrase VIII;Cals1;Car8;wdl
Protein Construction:	A DNA sequence encoding the mouse CA8 (P28651) (Met 1-Gln 291) was expressed, with a C-terminal polyhistidine tag. Predicted N terminal: Met 1
Species:	Mouse
Expression Host:	E. coli
Accession:	P28651
Molecular Weight:	34.5 kDa (predicted); 37 kDa (reducing conditions)

QC Testing

Biological Activity:	Measured by its esterase activity. The specific activity is >5 pmoles/min/ μ g.
Purity:	> 88 % as determined by SDS-PAGE
Endotoxin:	Please contact us for more information.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing 50 mM Tris, pH 8.0. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

The carbonic anhydrases (or carbonate dehydratases) are classified as metalloenzyme for its zinc ion prosthetic group and form a family of enzymes that catalyze the rapid interconversion of carbon dioxide and water to bicarbonate and protons, a reversible reaction that takes part in maintaining acid-base balance in blood and other tissues. The carbonic anhydrase (CA) family consists of at least 11 enzymatically active members and a few inactive homologous proteins. Carbonic anhydrase protein (CA) VIII, which is a member of the CA gene family, has

been shown to have no catalytic CA activity and its biological function is still unknown. Increased expression of CA-RP VIII was observed in 78% of colorectal carcinomas. It suggested that CA-RP VIII plays a role in the process of invasion in colorectal cancer.

Reference

Lehtonen J, et al. (2004) Characterization of CA XIII, a Novel Member of the Carbonic Anhydrase Isozyme Family. The Journal of Biological Chemistry. 279: 2719-27.

Lindskog S. (1997) Structure and mechanism of carbonic anhydrase. Pharmacology & Therapeutics. 74(1): 1-20.

Baird TT, et al. (1997) Catalysis and Inhibition of Human Carbonic Anhydrase IV. Biochemistry. 36 (9): 2669-78.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481