

ALDH4A1 Protein, Human, Recombinant (His & GST)

General Information

Synonyms:	P5CD;aldehyde dehydrogenase 4 family, member A1;ALDH4;P5CDh
Protein Construction:	A DNA sequence encoding the mature form of human ALDH4A1 (AAH07581.1) (Lys 25-Gln 563) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus. Predicted N terminal: Met
Species:	Human
Expression Host:	Baculovirus Insect Cells
Accession:	AAH07581.1
Molecular Weight:	86.8 kDa (predicted); 75 kDa (reducing conditions)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 80 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 20 mM Tris, 500 mM NaCl, pH 8.5, 10% glycerol. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

ALDH4A1 is a member of the aldehyde dehydrogenase family. Aldehyde dehydrogenase enzymes function in the metabolism of many molecules including certain fats (cholesterol and other fatty acids) and protein building blocks (amino acids). Additional aldehyde dehydrogenase enzymes detoxify external substances, such as alcohol and pollutants, and internal substances, such as toxins that are formed within cells. ALDH4A1 is expressed

abundantly in liver followed by skeletal muscle, kidney, heart, brain, placenta, lung and pancreas. It is a mitochondrial matrix NAD-dependent dehydrogenase which catalyzes the second step of the proline degradation pathway, converting pyrroline-5-carboxylate to glutamate. Defects in ALDH4A1 are the cause of hyperprolinemia type 2 (HP-2). HP-2 is characterized by the accumulation of delta-1-pyrroline-5-carboxylate (P5C) and proline. The disorder may be causally related to neurologic manifestations, including seizures and mental retardation.

Reference

Goodman SI, et al. (1974) Defective hydroxyproline metabolism in type II hyperprolinemia. *Biochemical medicine*. 10 (4): 329-36.

Maruyama K, et al. (1994) Oligo-capping: a simple method to replace the cap structure of eukaryotic mRNAs with oligoribonucleotides. *Gene*. 138 (1-2): 171-4.

Vasiliou V, et al. (2005) Analysis and update of the human aldehyde dehydrogenase (ALDH) gene family. *Hum Genomics*. 2 (2): 138-43.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481