

IFN-alpha 13/IFNA13 Protein, Cynomolgus, Recombinant (hFc)

General Information

Synonyms:	IFN- α 13;interferon, alpha 13;IFNA13;interferon, α 13
Protein Construction:	A DNA sequence encoding the cynomolgus IFNA13 (G7NFW4) (Met1-Glu190) was expressed with the Fc region of human IgG1 at the C-terminus. Predicted N terminal: Cys 25
Species:	Cynomolgus
Expression Host:	HEK293 Cells
Accession:	G7NFW4
Molecular Weight:	46.4 kDa (predicted); 47 kDa (reducing conditions)

QC Testing

Biological Activity:	Measured in antiviral assay using WISH human amnion cells infected with vesicular stomatitisvirus (VSV). The ED50 for this effect is typically 0.5-3 ng/mL.
Purity:	> 85 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Reference

Mahmutovic S, et al.. (2004) Significance of the interferon (IFN) in the therapy. *Bosn J Basic Med Sci.* 4(4): 42-4.

Wang, et al.. (2004) Fever of recombinant human interferon-alpha is mediated by opioid domain interaction with opioid receptor inducing prostaglandin E2. *J Neuroimmunol.* 156(1-2): 107-12.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481