

KEAP1 Protein, Human, Recombinant

General Information

Synonyms:	kelch like ECH associated protein 1; INRF2; KLHL19; KEAP-1
Protein Construction:	A DNA sequence encoding the human KEAP1 (Q14145) (Gln2-Cys624) was expressed and purified with two additional amino acids (Gly & Pro) at the N-terminus. Predicted N terminal: Gly
Species:	Human
Expression Host:	Baculovirus Insect Cells
Accession:	Q14145
Molecular Weight:	69.7 kDa (predicted); 64 kDa (reducing condition, due to glycosylation)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 20 mM Tris, 500 mM NaCl, 3 mM DTT, 10% glycerol, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Kelch-like ECH-associated protein 1, also known as a cytosolic inhibitor of Nrf2, Kelch-like protein 19, KEAP1, and INRF2, is a cytoplasm and nucleus protein that contains one BACK (BTB/Kelch associated) domain, one BTB (POZ) domain, and six Kelch repeats. KEAP1 / INRF2 is broadly expressed, with the highest levels in skeletal muscle. KEAP1 / INRF2 is a key regulator of the NRF2 transcription factor, which transactivates the antioxidant response

element (ARE) and upregulates numerous proteins involved in antioxidant defense. Under basal conditions, KEAP1 / INRF2 targets NRF2 for ubiquitination and proteolytic degradation and as such is responsible for the rapid turnover of NRF2. KEAP1 / INRF2 retains NFE2L2 / NRF2 in the cytosol. KEAP1 / INRF2 functions as a substrate adapter protein for the E3 ubiquitin ligase complex formed by CUL3 and RBX1. It targets NFE2L2 / NRF2 for ubiquitination and degradation by the proteasome, thus resulting in the suppression of its transcriptional activity and the repression of antioxidant response element-mediated detoxifying enzyme gene expression. KEAP1 / INRF2 may also retain BPTF in the cytosol. It targets PGAM5 for ubiquitination and degradation by the proteasome.

Reference

Strachan GD. et al., 2004, Biochemistry. 43: 12113-22.

Wang J, Chen B, Cheng C, et al. Timosaponin B II as a novel KEAP1-NRF2 inhibitor to alleviate alcoholic liver disease: Receptor structure-based virtual screening and biological evaluation. *Chemico-Biological Interactions*. 2025: 111390.

Zhang DD. et al., 2004, Mol Cell Biol. 24: 10941-53.

Zhang DD. et al., 2003, Mol Cell Biol. 23: 8137-51.

Lo S-C. et al., 2006, J Biol Chem. 281: 37893-903.

Buckley BJ. et al., 2008, Free Radic Biol Med. 44 (4) :692-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481