

## Endothelin B Receptor Protein, Human, Recombinant (mFc)

### General Information

Synonyms:	ET-B;endothelin receptor type B;ETB;ETBR;ETRB;ET-BR;HSCR;WS4A;HSCR2;ABCDS
Protein Construction:	A DNA sequence encoding the human EDNRB (NP_000106.1) (Met1-Lys101) was expressed with the Fc region of mouse IgG1 at the C-terminus. Predicted N terminal: Glu 27
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P24530-1
Molecular Weight:	34.4 kDa (predicted); 47.5 and 43.7 kDa (reducing conditions)

### QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

The hypermethylation of EDNRB gene was remarkably related to infiltration and metastasis of gastric cancer and may attribute to the tumor progression. EDNRB is a new candidate tumor suppressor gene which is often down-regulated or even silenced by promoter hypermethylation in various human cancers. Low EDNRB expression played a role in the progression of ACC tumors. The autosomal recessive mutation in EDNRB may underlie a part of WS1 with the current diagnostic criteria, and supported that Hirschsprung's disease is a multifactorial genetic

disease which requires additional factors.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481