

CK2 alpha/CSNK2A1 Protein, Mouse, Recombinant (His & GST)

General Information

| | |
|-----------------------|--|
| Synonyms: | casein kinase 2, alpha 1 polypeptide;casein kinase 2, α 1 polypeptide;CK2 α /CSNK2A1; Csnk2a1-rs4 |
| Protein Construction: | A DNA sequence encoding the mouse CSNK2A1 (Q60737) (Met1-Gln391) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus. Predicted N terminal: Met |
| Species: | Mouse |
| Expression Host: | Baculovirus Insect Cells |
| Accession: | Q60737 |
| Molecular Weight: | 72.9 kDa (predicted); 65 kDa (reducing conditions) |

QC Testing

| | |
|----------------------|--|
| Biological Activity: | The specific activity was determined to be 34 nmol/min/mg using casein as substrate. |
| Purity: | > 85 % as determined by SDS-PAGE |
| Endotoxin: | < 1.0 EU/ μ g of the protein as determined by the LAL method. |
| Formulation: | Supplied as sterile 20 mM Tris, 500 mM NaCl, pH 8.0, 10% gly. |

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store the product under sterile conditions at -20°C to -80°C . Samples are stable for up to 12 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

Proteins are shipped with blue ice.

Protein Background

Casein kinase II subunit alpha, also known as CK II alpha, CSNK2A1 and CK2A1, is a member of the protein kinase superfamily, Ser / Thr protein kinase family and CK2 subfamily. Casein kinase II (CSNK2A1) is a serine / threonine protein kinase that phosphorylates acidic proteins such as casein. This kinase is composed of an alpha, an alpha-prime, and two beta subunits. The alpha subunits contain the catalytic activity while the beta subunits undergo autophosphorylation. Casein kinase II (CSNK2A1) is a constitutively active, ubiquitously expressed serine / threonine protein kinase that is thought to have a regulatory function in cell proliferation, cell differentiation and apoptosis. CSNK2A1 functions as a tetrameric complex consisting of two regulatory beta-subunits and two catalytic units (alpha and alpha') in a homomeric or heteromeric conformation. Whilst the alpha- and alpha'-

subunits are catalytically identical, proteins that regulate CSNK2A1, such as cdc2 and Hsp90, preferentially bind to the alpha and not the alpha'-subunit. CSNK2A1 can phosphorylate a number of key intracellular signaling proteins implicated in tumor suppression (p53 and PTEN) and tumorigenesis (myc, jun, NF-kappaB). CSNK2A1 is also thought to influence Wnt signaling via beta-catenin phosphorylation and the PI 3-K signaling pathway via the phosphorylation of Akt.

Reference

- Schlpfer J, et al. (1997) A radiation hybrid framework map of bovine chromosome 13. *Chromosome Res.* 5(8): 511-9.
- Wirkner U, et al. (1994) The human gene (CSNK2A1) coding for the casein kinase II subunit alpha is located on chromosome 20 and contains tandemly arranged Alu repeats. *Genomics.* 19(2): 257-65.
- Wirkner U, et al. (1998) Genomic organization and promoter identification of the human protein kinase CK2 catalytic subunit alpha (CSNK2A1). *Genomics.* 48(1): 71-8.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481