

FSH Beta Protein, Human, Recombinant (hFc)

General Information

Synonyms:	follicle stimulating hormone, beta polypeptide; follicle stimulating hormone, β polypeptide; FSH β ; HH24
Protein Construction:	A DNA sequence encoding the human FSHB (NP_000501.1) (Met1-Glu129) was expressed with the Fc region of human IgG1 at the C-terminus. Predicted N terminal: Asn 19
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P01225
Molecular Weight:	39.2 kDa (predicted)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95 % as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/ μ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C . For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Follicle-stimulating hormone (FSH) is a gonadotrope-derived heterodimeric glycoprotein. FSH plays an essential role in processes involved in human reproduction, including spermatogenesis and the ovarian cycle. FSHB represents a conservative vertebrate gene with a unique function and it is located in a structurally stable gene-poor region. Polymorphisms in the follicle stimulating hormone beta subunit (FSHB) and follicle stimulating

A DRUG SCREENING EXPERT

hormone receptor (FSHR) genes might disturb normal spermatogenesis and affect male reproductive ability. The FSHB -211G>T genotype is a key determinant in the regulation of gonadotropins in different reproductive-endocrine pathophysiology.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481