

PLA2G2A Protein, Mouse, Recombinant (His)

General Information

Synonyms:	Pla2;sPLA2;sPla2-IIA;Mom1;EF;phospholipase A2, group IIA (platelets, synovial fluid)
Protein Construction:	A DNA sequence encoding the mouse PLA2G2A (NP_001076000.1) (Asn22-Cys146) was expressed with a polyhistidine tag at the C-terminus. Predicted N terminal: Asn 22
Species:	Mouse
Expression Host:	P. pastoris (Yeast)
Accession:	P31482
Molecular Weight:	15.3 kDa (predicted)

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95 % as determined by SDS-PAGE.
Endotoxin:	Please contact us for more information.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 50 mM NaAC, pH 4. 0. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Phospholipase A2, membrane associated, also known as Phosphatidylcholine 2-acylhydrolase 2A, Group IIA phospholipase A2, Non-pancreatic secretory phospholipase A2 and PLA2G2A, is a peripheral membrane protein that belongs to the phospholipase A2 family. PLA2G2A is found in many cells and also extracellularly. The membrane-bound and secreted forms of PLA2G2A are identical. PLA2G2A has been proposed to play a role in anti-bacterial defense, inflammation and eicosanoid generation, in clearance of apoptotic cells, and in the Wnt

signaling pathway. PLA2G2A is thought to participate in the regulation of the phospholipid metabolism in biomembranes including eicosanoid biosynthesis. PLA2G2A catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides. PLA2G2A might be a factor in human colorectal tumorigenesis.

Reference

- Praml,C. et al., 1998, Oncogene. 17 (15):2009-12.
Fijneman,R.J. et al., 2008,Front Biosci 13 :4144-74.
Fijneman,R.J. et al., 2009, Cell Oncol 31 (5):345-56.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481