

## ICA Protein, Mouse, Recombinant (His)

### General Information

Protein Construction:	A DNA sequence encoding the mouse ICA (NP_082194.1) (Met1-Tyr700) was expressed with a polyhistidine tag at the C-terminus. Predicted N terminal: Leu 20
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	Q9DBD0
Molecular Weight:	76.2 kDa (predicted)

### QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95 % as determined by SDS-PAGE
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 $\mu$ m filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

INHCA, official full name is inhibitor of carbonic anhydrase, which is a member of the transferrin gene family. Unlike many members of that gene family it does not play a significant role in iron homeostasis but rather is an inhibitor of carbonic anhydrase. Carbonic anhydrase inhibitors are a class of pharmaceuticals that suppress the activity of carbonic anhydrase.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481