

TRIM21 Protein, Human, Recombinant (His)

General Information

| | |
|-----------------------|--|
| Synonyms: | Ro/SSA;SSA;RO52;RNF81;tripartite motif containing 21;SSA1 |
| Protein Construction: | A DNA sequence encoding the human TRIM21 (NP_003132.2) (Met1-Tyr475) was expressed with a polyhistidine tag at the N-terminus. Predicted N terminal: His |
| Species: | Human |
| Expression Host: | Baculovirus Insect Cells |
| Accession: | P19474-1 |
| Molecular Weight: | 56.5 kDa (predicted) |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first. |
| Purity: | > 90 % as determined by SDS-PAGE. |
| Endotoxin: | < 1.0 EU/µg of the protein as determined by the LAL method. |
| Formulation: | Lyophilized from a solution filtered through a 0.22 µm filter, containing 20 mM Tris, pH 8.0, 500 mM NaCl, 10% glycerol, 1 mM TCEP. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization. |

Preparation and Storage

| | |
|----------------------|--|
| Reconstitution: | Reconstituted with sterile deionized water to 0.25 mg/mL. Reconstitution conditions may vary depending on the lot. |
| Stability & Storage: | It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small> |
| Shipping: | In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice. |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481