

## CD40 Protein, Human, Recombinant (His & Avi), Biotinylated

### General Information

Synonyms:	Bp50;p50;CDW40;CD40 molecule, TNF receptor superfamily member 5;TNFRSF5
Protein Construction:	A DNA sequence encoding the human CD40 (NP_001241.1) (Met1-Arg193) was expressed with a c-terminal polyhistidine tagged AVI tag at the C-terminus. The expressed protein was biotinylated in vivo by the Biotin-Protein ligase (BirA enzyme) which is co-expressed. Predicted N terminal: Glu 21
Species:	Human
Expression Host:	HEK293 Cells
Accession:	P25942
Molecular Weight:	22.4 kDa (predicted)

### QC Testing

Biological Activity:	Immobilized Human CD40 Protein (ECD, His & AVI Tag), Biotinylated at 2 µg/mL (100 µL/well) can bind Human CD40 Ligand Protein (ECD, hFc Tag), the EC50 is 6-40 ng/mL.
Purity:	> 95 % as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

Reconstitution:	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

CD40, also known as TNFRSF5, is a member of the TNF receptor superfamily which are single transmembrane-spanning glycoproteins. CD40 protein plays an essential role in mediating a broad variety of immune and inflammatory responses including T cell-dependent immunoglobulin class switching, memory B cell development,

and germinal center formation. CD40 protein is expressed in B cells, dendritic cells, macrophages, endothelial cells, and several tumor cell lines. Defects in CD40 result in hyper-IgM immunodeficiency type 3 (HIGM3). In addition, CD40/CD40L interaction is found to be necessary for amyloid-beta-induced microglial activation, and thus is thought to be an early event in Alzheimer disease pathogenesis. Cancer Immunotherapy Co-stimulatory Immune Checkpoint Targets Immune Checkpoint Immune Checkpoint Detection: Antibodies Immune Checkpoint Detection: ELISA Antibodies Immune Checkpoint Detection: FCM Antibodies Immune Checkpoint Detection: ICC Antibodies Immune Checkpoint Detection: IP Antibodies Immune Checkpoint Detection: WB Antibodies Immune Checkpoint Proteins Immune Checkpoint Targets Immunotherapy Targeted Therapy

### Reference

- van Kooten C, et al. (2000). CD40-CD40 ligand. *J Leukoc Biol.* 67 (1): 2-17.
- Bhushan A, et al. (2002). CD40:CD40L interactions in X-linked and non-X-linked hyper-IgM syndromes. *Immunol Res.* 24 (3): 311-24.
- Chatzigeorgiou A, et al. (2009) CD40/CD40L signaling and its implication in health and disease. *Biofactors.* 35(6): 474-83.
- Li R, et al. (2009) Expression of CD40 and CD40L in Gastric Cancer Tissue and Its Clinical Significance. *Int J Mol Sci.* 10 (9): 3900-17.
- Lievens D, et al. (2009) The multi-functionality of CD40L and its receptor CD40 in atherosclerosis. *Thromb Haemost.* 102(2): 206-14.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481