

Influenza A H3N2 (A/Hong Kong/2671/2019) Nucleoprotein/NP Protein (His)

General Information

Synonyms:	NP Protein
Protein Construction:	A DNA sequence encoding the Influenza A virus (A/Hong Kong/2671/2019 (H3N2)) nucleoprotein (translated amino acids of EPI1698482) (Met1-Asn498) was expressed with a polyhistidine tag at the C-terminus. Predicted N terminal: Met 1
Species:	H3N2
Expression Host:	Baculovirus Insect Cells
Molecular Weight:	57.3 kDa (predicted)

QC Testing

Biological Activity:	Activity testing is in progress. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95 % as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/μg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 μm filter, containing 20 mM Tris, 500 mM NaCl, 10% glycerol, pH 8.0. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

Preparation and Storage

Reconstitution:	Reconstituted with 20 mM Tris, 500 mM NaCl, pH 8.0 to 0.25 mg/mL. Reconstitution conditions may vary depending on the lot.
Stability & Storage:	It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Protein Background

Influenza viral nucleoprotein (NP) plays a critical role in virus replication and host adaptation. Nucleoprotein (NP) is an essential multifunctional protein that encapsidates the viral genome and functions as an adapter between the virus and the host cell machinery. NPs contain two nuclear localization signals (NLSs): a well-studied monopartite NLS1 and a less-characterized NLS2, thought to be bipartite. The nucleocapsid is a complex of the viral nucleoprotein, RNA, and several other viral proteins. The nucleoprotein forms large, RNA-bound, helical

filaments and acts as a scaffold for additional viral proteins.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481