

## BAFFR/TNFRSF13C Protein, Human, Recombinant (hFc &amp; Avi), Biotinylated

## General Information

|                       |   |
|-----------------------|---|
| Synonyms:             | BAFFR;tumor necrosis factor receptor superfamily, member 13C;BROMIX;prolixin;BAFF-R;CD268;CVID4   |
| Protein Construction: | A DNA sequence encoding the human TNFRSF13C (NP_443177.1)(Arg2-Ala71) was expressed with a c-terminal AVI tagged Fc region of human IgG1 at the C-terminus (Fc-AVI). The expressed protein was biotinylated in vivo by the Biotin-Protein ligase (BirA enzyme) which is co-expressed. Predicted N terminal: Ser |
| Species:              | Human   |
| Expression Host:      | HEK293 Cells  |
| Accession:            | Q5H8V1  |
| Molecular Weight:     | 35.79 kDa (predicted); 51.85 and 44.65 kDa (reducing conditions)  |

## QC Testing

|                      |  |
|----------------------|--|
| Biological Activity: | Immobilized Human BAFF/BlyS hFc at 2 µg/mL (100 µL/well) can bind Human BAFFR/TNFRSF13C hFc & AVI, Biotinylated, the EC50 of Human BAFFR/TNFRSF13C hFc & AVI, Biotinylated is 40-200 ng/mL.  |
| Purity:              | > 95 % as determined by SDS-PAGE.  |
| Endotoxin:           | < 1.0 EU/µg of the protein as determined by the LAL method.  |
| Formulation:         | Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization. |

## Preparation and Storage

|                      |  |
|----------------------|--|
| Reconstitution:      | A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.  |
| Stability & Storage: | It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.<br><small>Actual storage temperature shall be subject to the COA.</small> |
| Shipping:            | In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.   |

## Protein Background

Tumor necrosis factor receptor superfamily, member 13C (TNFRSF13C) also known as B-cell-activating factor

receptor (BAFFR) and CD268 antigen, is a member of the tumor necrosis factor receptor superfamily. A tumor necrosis factor receptor (TNFR), or death receptor, is a trimeric cytokine receptor that binds tumor necrosis factors (TNF). The receptor cooperates with an adaptor protein which is important in determining the outcome of the response. Members of the TNF receptor superfamily (TNFRSF) have crucial roles in both innate and adaptive immunity and in cellular apoptosis process. Apoptosis is a cell suicide mechanism that enables metazoans to control cell number in tissues and to eliminate individual cells that threaten the animal's survival. Certain cells have unique sensors, termed death receptors or tumour necrosis factor (TNFR), on their surface. Tumour necrosis factors (TNFR) detect the presence of extracellular death signals and, in response, they rapidly ignite the cell's intrinsic apoptosis machinery. It has been proposed that abnormally high levels of BAFFR/TNFRSF13C (CD268) may contribute to the pathogenesis of autoimmune diseases by enhancing the survival of autoreactive B cells.

### Reference

- Ashkenazi A, et al. (1998) Death receptors: signaling and modulation. *Science*. 281(5381): 1305-8.
- Losi CG, et al. (2005) Mutational analysis of human BAFF receptor TNFRSF13C (BAFF-R) in patients with common variable immunodeficiency. *J Clin Immunol*. 25(5): 496-502.
- Hentges KE, et al. (2002) Tnfrsf13c (Baffr) is mis-expressed in tumors with murine leukemia virus insertions at Lvis22. *Genomics*. 80(2): 204-12.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481