

## BTN1A1 Protein, Human, Recombinant (His & Avi), Biotinylated

### General Information

Synonyms:	BTN;BT;BTN1
Protein Construction:	Ala27-Arg242
Species:	Human
Expression Host:	HEK293 Cells
Accession:	Q13410
Molecular Weight:	26.9 kDa (Predicted); 30-40 kDa (Reducing conditions due to glycosylation)

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it.
Purity:	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC
Endotoxin:	< 1.0 EU/ $\mu$ g of the protein as determined by the LAL method.
Formulation:	Lyophilized from 0.22 $\mu$ m filtered solution in PBS (pH 7.4). Normally 8% trehalose is added as protectant before lyophilization.

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100  $\mu$ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

Butyrophilin 1A1 (BTN1A1) is one of the membrane proteins that surrounds LD, BTN1A1 plays an important role in regulating LD synthesis via a mechanism involving membrane phospholipid composition.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481