

## Monkeypox Virus (MPXV) Envelope protein E8L Protein (His)

## General Information

Protein Construction:	A DNA sequence encoding the Monkeypox Virus (MPXV) E8L Protein (QNP13700.1) (Met1-Thr275) was expressed with a polyhistidine tag at the C-terminus. Predicted N terminal: Met 1
Species:	MPXV
Expression Host:	Baculovirus Insect Cells
Accession:	QNP13700.1
Molecular Weight:	35.75 kDa (predicted); 67.64 kDa and 37.35 kDa (reducing conditions)

## QC Testing

Biological Activity:	Immobilized Monkeypox Virus (MPXV) Envelope protein E8L Antibody, Mouse MAb at 2 µg/mL (100 µL/well) can bind Recombinant Monkeypox Virus (MPXV) Envelope protein E8L Protein (His Tag) (Cat#TMPY-06562), the EC50 is 5-20 ng/mL.
Purity:	> 90 % as determined by SDS-PAGE.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing 20 mM Tris, 300 mM NaCl, 10% glycerol, pH 8.5. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

## Preparation and Storage

## Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

## Stability &amp; Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481