

## B7-1 Protein, Mouse, Recombinant (His & Avi), Biotinylated

### General Information

Synonyms:	MIC17;Ly-53;Cd28l;B71;Ly53;TSA1;CD80 molecule
Protein Construction:	A DNA sequence encoding the Mouse CD80 (Q00609-1) (Met1-Lys245) was expressed with a C-terminal polyhistidine tag followed by an AVI tag. The expressed protein was biotinylated in vivo by the Biotin-Protein ligase (BirA enzyme) which is co-expressed. Predicted N terminal: Val 38
Species:	Mouse
Expression Host:	HEK293 Cells
Accession:	Q00609-1
Molecular Weight:	26.94 kDa (predicted); 56.66 kDa (reducing conditions)

### QC Testing

Biological Activity:	Immobilized Recombinant Mouse CTLA4 / CD152 Protein (ECD, Fc Tag) at 2µg/mL (100µL/well) can bind Recombinant Mouse B7-1 Protein (ECD, His & AVI Tag), Biotinylated , the EC50 is 12-32 ng/mL.
Purity:	≥ 95 % as determined by SDS-PAGE. ≥ 90 % as determined by SEC-HPLC.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

**Reconstitution:**  
A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

The B-lymphocyte activation antigen B7-1 (referred to as B7), also known as CD80, is a member of cell surface immunoglobulin superfamily and is expressed on the surface of antigen-presenting cells including activated B

cells, macrophages and dendritic cells. As costimulatory ligands, B7-1 which exists predominantly as dimer and the related protein B7-2, interact with the costimulatory receptors CD28 and cytotoxic T lymphocyte-associated antigen 4 (CTLA-4) expressed on T cells, and thus constitute one of the dominant pathways that regulate T cell activation and tolerance, cytokine production, and the generation of CTL. The B7/CD28/CTLA4 pathway has the ability to both positively and negatively regulate immune responses. CD80 is thus regarded as promising therapeutic targets for autoimmune diseases and various carcinomas. Cancer Immunotherapy Co-inhibitory Immune Checkpoint Targets Immune Checkpoint Immune Checkpoint Detection: Antibodies Immune Checkpoint Detection: ELISA Antibodies Immune Checkpoint Detection: FCM Antibodies Immune Checkpoint Detection: IHC Antibodies Immune Checkpoint Detection: WB Antibodies Immune Checkpoint Proteins Immune Checkpoint Targets Immunotherapy Targeted Therapy

### Reference

- Greenfield EA, et al. (1998) CD28/B7 costimulation: a review. *Crit Rev Immunol.* 18(5): 389-418.
- Zang X, et al. (2007) The B7 family and cancer therapy: costimulation and coinhibition. *Clin Cancer Res.* 13(18 Pt 1): 5271-9.
- Mir MA, et al. (2008) Signaling through CD80: an approach for treating lymphomas. *Expert Opin Ther Targets.* 12(8): 969-79.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481