

## B7-1 Protein, Rat, Recombinant (His)

### General Information

Synonyms:	CD80 molecule
Protein Construction:	A DNA sequence encoding the Rat CD80 (AAB60503.1) (Met1-Gln248) was expressed with a polyhistidine tag at the C-terminus. Predicted N terminal: Ile 39
Species:	Rat
Expression Host:	HEK293 Cells
Accession:	AAB60503.1
Molecular Weight:	25.38 kDa (predicted); 42.0 kDa (reducing conditions)

### QC Testing

Biological Activity:	Immobilized Recombinant Rat CTLA-4 / CD152 Protein (ECD, Fc Tag) at 2 µg/mL (100 µL/well) can bind Recombinant Rat B7-1 Protein (ECD, His Tag), the EC50 is 1.2-4.5 ng/mL.
Purity:	≥ 95 % as determined by SDS-PAGE. ≥ 95 % as determined by SEC-HPLC.
Endotoxin:	< 1.0 EU/µg of the protein as determined by the LAL method.
Formulation:	Lyophilized from a solution filtered through a 0.22 µm filter, containing PBS, pH 7.4. Typically, a mixture containing 5% to 8% trehalose, mannitol, and 0.01% Tween 80 is incorporated as a protective agent before lyophilization.

### Preparation and Storage

#### Reconstitution:

A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.

#### Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

### Protein Background

The B-lymphocyte activation antigen B7-1 (referred to as B7), also known as CD80, is a member of cell surface immunoglobulin superfamily and is expressed on the surface of antigen-presenting cells including activated B cells, macrophages and dendritic cells. As costimulatory ligands, B7-1 which exists predominantly as dimer and the related protein B7-2, interact with the costimulatory receptors CD28 and cytotoxic T lymphocyte-associated antigen 4 (CTLA-4) expressed on T cells, and thus constitute one of the dominant pathways that regulate T cell

activation and tolerance, cytokine production, and the generation of CTL. The B7/CD28/CTLA4 pathway has the ability to both positively and negatively regulate immune responses. CD80 is thus regarded as promising therapeutic targets for autoimmune diseases and various carcinomas. Cancer Immunotherapy Co-inhibitory Immune Checkpoint Targets Immune Checkpoint Immune Checkpoint Detection: Antibodies Immune Checkpoint Detection: ELISA Antibodies Immune Checkpoint Detection: FCM Antibodies Immune Checkpoint Detection: IHC Antibodies Immune Checkpoint Detection: WB Antibodies Immune Checkpoint Proteins Immune Checkpoint Targets Immunotherapy Targeted Therapy

### Reference

Greenfield EA, et al. (1998) CD28/B7 costimulation: a review. Crit Rev Immunol. 18(5): 389-418.

Zang X, et al. (2007) The B7 family and cancer therapy: costimulation and coinhibition. Clin Cancer Res. 13(18 Pt 1): 5271-9.

Mir MA, et al. (2008) Signaling through CD80: an approach for treating lymphomas. Expert Opin Ther Targets. 12(8): 969-79.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481