

Taxcultine

Chemical Properties


CAS No. : 153415-46-4

Formula: C44H53NO14

Molecular Weight: 819.89

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Taxcultine can actively promote tubulin assembly.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2197 mL	6.0984 mL	12.1968 mL
5 mM	0.2439 mL	1.2197 mL	2.4394 mL
10 mM	0.122 mL	0.6098 mL	1.2197 mL
50 mM	0.0244 mL	0.122 mL	0.2439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481