

Compound N008-0008

Chemical Properties


CAS No. :

Formula: C45H74O19

Molecular Weight: 919.068

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N008-0008 is a useful organic compound for research related to life sciences and the catalog number is T131369.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0881 mL	5.4403 mL	10.8806 mL
5 mM	0.2176 mL	1.0881 mL	2.1761 mL
10 mM	0.1088 mL	0.544 mL	1.0881 mL
50 mM	0.0218 mL	0.1088 mL	0.2176 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481