

Compound N014-0003

Chemical Properties


CAS No. :

Formula: C₃₆H₅₀O₁₅

Molecular Weight: 722.781

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N014-0003 is a useful organic compound for research related to life sciences and the catalog number is T131358.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3835 mL	6.9177 mL	13.8355 mL
5 mM	0.2767 mL	1.3835 mL	2.7671 mL
10 mM	0.1384 mL	0.6918 mL	1.3835 mL
50 mM	0.0277 mL	0.1384 mL	0.2767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481