

1-(b-D-Xylofuranosyl)-N6-isopentenyl adenine

Chemical Properties

CAS No. :


Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-(b-D-Xylofuranosyl)-N6-isopentenyl adenine is a Nucleoside Derivative - Xylo-nucleoside, 6-Modified purine nucleoside.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481