

Compound N001-0004

## Chemical Properties


CAS No. :

Formula: C<sub>29</sub>H<sub>44</sub>O<sub>12</sub>

Molecular Weight: 584.659

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N001-0004 is a useful organic compound for research related to life sciences and the catalog number is T131388.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7104 mL	8.552 mL	17.104 mL
5 mM	0.3421 mL	1.7104 mL	3.4208 mL
10 mM	0.171 mL	0.8552 mL	1.7104 mL
50 mM	0.0342 mL	0.171 mL	0.3421 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481