


Withaphysalin E

Chemical Properties

CAS No. : 118985-24-3
 Formula: C₂₈H₃₄O₇
 Molecular Weight: 482.57
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Withaphysalin E is a natural product used in life sciences research [catalog number TN5320, CAS number 118985-24-3].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0722 mL	10.3612 mL	20.7224 mL
5 mM	0.4144 mL	2.0722 mL	4.1445 mL
10 mM	0.2072 mL	1.0361 mL	2.0722 mL
50 mM	0.0414 mL	0.2072 mL	0.4144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481