


## Thieno[3,2-d]pyrimidin-4-amine, N,N-diethyl-

## Chemical Properties

CAS No. :	16234-52-9
Formula:	C10H13N3S
Molecular Weight:	207.3
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Thieno[3,2-d]pyrimidin-4-amine, N,N-diethyl-, with CAS No. 16234-52-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Thieno[3,2-d]pyrimidin-4-amine, N,N-diethyl- provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8239 mL	24.1196 mL	48.2393 mL
5 mM	0.9648 mL	4.8239 mL	9.6479 mL
10 mM	0.4824 mL	2.412 mL	4.8239 mL
50 mM	0.0965 mL	0.4824 mL	0.9648 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481