

4,6,8-trihydroxy-3-methyl-3,4-dihydro-1H-2-benzopyran-1-one

Chemical Properties


CAS No. :

Formula: C₁₀H₁₀O₅

Molecular Weight: 210.19

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

4,6,8-trihydroxy-3-methyl-3,4-dihydro-1H-2-benzopyran-1-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7576 mL	23.788 mL	47.576 mL
5 mM	0.9515 mL	4.7576 mL	9.5152 mL
10 mM	0.4758 mL	2.3788 mL	4.7576 mL
50 mM	0.0952 mL	0.4758 mL	0.9515 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481