

## 2-Benzylsuccinic acid


## Chemical Properties

CAS No. : 884-33-3

Formula: C<sub>11</sub>H<sub>12</sub>O<sub>4</sub>

Molecular Weight: 208.21

Storage: Keep away from moisture, Keep away from direct sunlight  
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	2-Benzylsuccinic acid (DL-Benzylsuccinic acid) is a metabolite produced by bacteria as they process xenobiotic compounds.
Targets(IC50)	Carboxypeptidase

## Solubility Information

Solubility	DMSO: 41 mg/mL (196.92 mM), Sonication is recommended. ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (19.21 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8028 mL	24.0142 mL	48.0284 mL
5 mM	0.9606 mL	4.8028 mL	9.6057 mL
10 mM	0.4803 mL	2.4014 mL	4.8028 mL
50 mM	0.0961 mL	0.4803 mL	0.9606 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481